

SARS-CoV-2 Antigen Test Kit (Colloidal Gold)

The COVID-19 is spread by human-to-human transmission and infection has been estimated to have a mean incubation period of 6.4 days and a basic reproduction number is 2.24-3.58. Many countries in the world have recognized the value of antigen detection in the early diagnosis of COVID-19 infection. The SARS-CoV-2 Antigen Test Kit is a gold immuno-chromatographic assay that is intended for the qualitative detection of specific antigen of SARS-CoV-2 in nasal swab, nasopharyngeal swab and saliva specimens of human. It can be used for self test at home or for professional test by a medical staff.

Key Feature



Cost-effective



No equipment required, simple operation



Rapid on-site test




Simple interpretation
of results




Fast results
within 15 minutes

Test Procedure



- 1 Open the aluminum foil bag, put the test cassette on a clean, horizontal bench. Place the swab into the sample tube that has been pre-filled with 0.5 ml sample buffer, rotate the swab for about 10 seconds, and press the swab applicator against the tube wall to release the antigen in the swab.
- 2 Roll the swab applicator against the inside of the sample tube as you remove it. Dispose of the used swab in a biohazard waste following the local government regulations.
- 3 Install the dropper cap onto the sample tube.
- 4 Add Two Drops of the extraction solution into the sample well and start the timer. Wait for 5~15 minutes and observe the results. Results are invalid if it appeared exceed 15 minutes.

Interpretation of results



Ordering Information

Item	Component	Specification Qty.	
		1 Test/kit R041-01	20 Tests/kit R041-02
1	Test Cassette	1	20
2	Sample Tube, with 0.5 ml sample buffer	1	20
3	Single packaged swab	1	20
4	Instruction for use	1	1

SHENZHEN UNI-MEDICA TECHNOLOGY CO., LTD

Address: Block 6th, Liuxian Culture Park, Nanshan District,
Shenzhen, Guangdong Province, P. R. China

Telephone: 0755-86502782

E-mail: marketing@uni-medica.com

SARS-CoV-2 Antigen Test Kit (Colloidal Gold)

INSTRUCTION FOR USE

INTENDED USE

The SARS-CoV-2 Antigen Test Kit is a gold immuno-chromatographic assay (GICA) that is intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal (NP) swab, nasal (NS) swab, and saliva specimens directly from individuals who are suspected of COVID-19 by their healthcare provider, with or without signs and symptoms of COVID-19.

The SARS-CoV-2 Antigen Test Kit does not differentiate between SARS-CoV and SARS-CoV-2.

Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. The antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of the disease.

If necessary, negative results should be treated as presumptive and confirmed with a molecular assay for patient management. Negative results do not rule out COVID-19 and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19.

The SARS-CoV-2 Antigen Test Kit is intended for use by trained clinical laboratory personnel and individuals trained in point of care settings.

GENERAL INFORMATION

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue, and dry cough. Nasal congestion, runny nose, sore throat, myalgia, and diarrhea are found in a few cases. The median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12 days of infection. The symptoms of COVID-19 are similar to other viral respiratory diseases including fever, cough, and shortness of breath.

PRINCIPLE OF THE TEST

The SARS-CoV-2 Antigen Test Kit is a rapid lateral flow immuno-chromatographic sandwich assay to directly detect nucleocapsid protein of SARS-CoV-2 in the nasopharyngeal swab, nasal swab, and saliva specimens and diagnosis of SARS-CoV-2 infection.

The patient sample is placed in the Sample Tube, during which time the virus particles in the sample are disrupted, exposing internal viral nucleoproteins. After disruption, the sample is added to the Test Cassette sample well. And the sample migrates through a test strip, if the SARS-CoV-2 virus antigen is present, a red color line will be shown on the T line. If the SARS-CoV-2 viral antigen is absent, there is not a red line that will be shown on the T line, however, a red line will be always shown on the C line indicating that the reaction system properly happens.

REAGENTS AND MATERIALS PROVIDED

Item	Components	Specification/Qty.	
		1 Test/kit	20 Tests/kit
1	Test Cassette individually foil pouched with a desiccant	1	20
2	Sample Tube, with 0.5 mL sample buffer.	1	20
3	Single packaged NP/NS swab	1	20
4	Instruction for use	1	1

*** Components will be included when customers demand.**

Materials needed but not provided:

- Timer or watch.
- Vortex
- Saliva collection device/cup/bag
- 1.0/0.3-mL transfer pipette

PRECAUTIONS

- For in vitro diagnostic use only.
- Please read this manual carefully before using this test kit. And follow the testing procedures strictly described in the manual. Otherwise, it will lead to incorrect results.
- Do not use expired reagents.
- Do not re-use the test kit.
- All swab samples, used reagents, test cards, and other materials used during testing are considered to be infectious, and personal protection should be done during the experiment.
- Use of Nitrile, Latex (or equivalent) gloves is recommended when handling patient samples. Wear suitable protective clothing and eye/face protection when handling the contents of this kit.
- Sample handling and waste disposal must comply with relevant regulations. Wash hands thoroughly after handling.
- Avoid using visually bloody or overly viscous samples for testing.
- Do not use components from different batch lots.
- The sample tube contains a salt solution. If the solution contacts the skin or eye, flush with copious amounts of water.
- Sample collection and handling procedures require specific training and guidance.

STORAGE AND STABILITY

- The test device is sensitive to humidity as well as heat.
- Store kit components at 2-30°C, out of direct sunlight. Kit components are stable until the expiration date printed on the outer box.
- After unsealing the aluminum foil bag, the test cassette should be used as soon as possible within **Two Hours**.
- Do not freeze.

SPECIMEN COLLECTION AND PREPARATION

Test specimens immediately after collection for optimal test performance. Inadequate specimen collection or improper sample handling/storage/transport may yield erroneous results. Refer to the CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19)

<https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>

- Nasal swab:

To collect a nasal swab sample, carefully insert the swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation, rotate the swab against the wall of the nostril into the patient's nostril to the nasal palate, and then slowly remove it while wiping.

- Nasopharyngeal swab:

To collect a nasopharyngeal swab sample, carefully insert the swab into the nostril that presents the most secretion under visual inspection. Keep the swab near the septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times then remove it from the nasopharynx.

- Saliva

Collect the saliva specimen 1 mL using a clean collecting device/cup/bag, then take 0.5 mL saliva sample into the sample tube.

SAMPLE TRANSPORT AND STORAGE

Samples should be tested as soon as possible after collection. The nasal or nasopharyngeal swabs are stable for up to 24 hours at room temperature or 2-8°C.

TEST PROCEDURE

Please read the instructions for use carefully before testing, and complete the test in strict accordance with the directions of the manual, otherwise reliable results cannot be guaranteed.

- Open the aluminum foil bag, put the test cassette on a clean, horizontal bench.
- Bring the samples to room temperature prior to assay in case of the samples were stored at 2-8°C.


Swab Test Procedure (Nasal/Nasopharyngeal):

- Place the swab into the sample tube that has been pre-filled with 0.5 mL sample buffer, rotate the swab for about 10 seconds, and press the swab applicator against the tube wall to release the antigen in the swab.
- Roll the swab applicator against the inside of the sample tube as you remove it. Dispose of the used swab in a biohazard waste following the local government regulations.
- Install the dropper cap onto the sample tube, add **Two Drops** of the extraction solution into the sample well and start the timer.

Saliva Test Procedure:

- Transfer 0.3 mL saliva specimen into the sample tube, vortex to extract the viral antigen in the specimen.
- Install the dropper cap onto the sample tube, add **Two Drops** of the extraction solution into the sample well and start the timer.

- Read the results within 15 minutes. And the results are invalid after 15 minutes.



INTERPRETATION OF TEST RESULTS

Negative result	Positive result	Invalid result	

1. Positive: A red line appears on the test line (T) and the control line (C).

NOTE: A positive result does not rule out co-infections with other pathogens.

2. Negative:

Only the control line (C) appears, and no red line appears on the test line (T).

NOTE: A negative result does not exclude infection.

3. Invalid:

There is no red line at the position of the control line (C). Regardless of whether the TEST line (T) is displayed, it is an invalid result and the sample should be tested again.

PERFORMANCES

LoD:

The LoD of the test kit is 10 pg/mL for detection with recombinant SARS-CoV-2 nucleocapsid protein and is 30 TCID₅₀ with inactive viral culture.

The Negative Percent Agreement (NPA):

The NPA of the test kit should be 20/20 (-/-) using an internal negative reference panel.

The Positive Percent Agreement (PPA):

The PPA of the test kit should be 8/8 (+/+) using an internal positive reference panel.

Repeatability:

The repeatability of the test kit should be tested using the same batch number, and all of the test results should be positive, and the T lines have even intensities.

ANALYTICAL PERFORMANCE

Limit of Detection (LoD)

The limit of detection (LoD) of the SARS-CoV-2 Antigen Test Kit was determined by evaluating different concentrations of inactivated SARS-CoV-2 virus from viral cell culture and of recombinant nucleocapsid protein (rNp). Tested nasal swab samples were prepared by absorbing 20 uL of each virus or rNp dilution onto the swab. The swab samples were tested according to the test procedure.

The LoD was determined as the lowest virus concentration that was detected $\geq 95\%$ of the time (i.e., the concentration at which at least 20 out of 21 replicates tested positive).

The LoD of the test kit in the nasal swab was confirmed as 30

TCID₅₀/swab and 10 pg/mL of rNp.

LoD Test Results

Concentration	Number: Positive/Total	Detected Percentage
30 TCID ₅₀ /swab	20/21	95%
10 pg/mL	20/21	95%

High Dose Hook Effect

No high dose hook effect was observed when tested with up to a concentration of 3.0x 10⁵ TCID₅₀/mL of inactivated SARS-CoV-2 virus with the test kit.

Cross-Reactivity

Cross-reactivity of SARS-CoV-2 Antigen Test Kit was evaluated by testing normal respiratory tract pathogenic microorganisms (bacteria, viruses, yeast) and a pooled human nasal wash that may be present in the nasal cavity. Each of the bacteria, viruses, and yeast were tested in duplicate in the absence or presence of heat-inactivated SARS-CoV-2 virus (1.5 x 10² TCID₅₀/swab). No cross-reactivity or interference came out of these tests when tested at the concentration listed below.

Potential Cross-Reactant		Test Concentration
Virus	Respiratory Syncytial Virus A	1.0 x 10 ⁵ PFU/mL
	Measles Virus	1.0 x 10 ⁵ TCID ₅₀ /mL
	Mumps Virus	1.0 x 10 ⁵ TCID ₅₀ /mL
	Adenovirus Type 3	1.0 x 10 ⁵ TCID ₅₀ /mL
	Parainfluenza Type 2	1.0 x 10 ⁵ TCID ₅₀ /mL
	Partial Pulmonary Virus	1.0 x 10 ⁵ TCID ₅₀ /mL
	Human coronavirus OC43	1.0 x 10 ⁵ TCID ₅₀ /mL
	Human coronavirus 229E	1.0 x 10 ⁵ TCID ₅₀ /mL
	Human coronavirus NL63	1.0 x 10 ⁵ TCID ₅₀ /mL
	MERS coronavirus	1.0 x 10 ⁵ TCID ₅₀ /mL
	Influenza B (Victoria Strain)	1.0 x 10 ⁵ TCID ₅₀ /mL
	Influenza B (Y Strain)	1.0 x 10 ⁵ TCID ₅₀ /mL
	Influenza A (H1N1,2009)	1.0 x 10 ⁵ TCID ₅₀ /mL
	Influenza A (H3N2)	1.0 x 10 ⁵ TCID ₅₀ /mL
	Avian Influenza Virus H7N9	1.0 x 10 ⁵ TCID ₅₀ /mL
	Avian Influenza Virus H5N1	1.0 x 10 ⁵ TCID ₅₀ /mL
Bacteria	Epstein Barr Virus	1.0 x 10 ⁵ TCID ₅₀ /mL
	Enterovirus CA16	1.0 x 10 ⁵ TCID ₅₀ /mL
	Rhinovirus	1.0 x 10 ⁵ PFU/mL
	Staphylococcus aureus	1.0 x 10 ⁶ CFU/mL
	Streptococcus pneumoniae	1.0 x 10 ⁶ CFU/mL
	Mycoplasma pneumoniae	1.0 x 10 ⁶ IFU/mL
	Parapertussis	1.0 x 10 ⁶ cells/mL
	Chlamydia pneumoniae	1.0 x 10 ⁶ IFU/mL
	Pooled human nasal wash	N/A

Yeast	<i>Candida albicans</i>	1.0 x 10 ⁶ cells/mL
-------	-------------------------	--------------------------------

CLINICAL PERFORMANCES

The clinical performance characteristics of the SARS-CoV-2 Antigen Test Kit were evaluated and a CE marked real-time Polymerase Chain Reaction (RT-PCR) assay for the detection of SARS-CoV-2 was used as the comparison method for this study.

SARS-CoV-2 Antigen Test Kit Performance against the RT-PCR Method

SARS-CoV-2 Antigen Test Kit	Comparison Method (RT-PCR)		
	Positive	Negative	Total
Positive	104	2	106
Negative	4	197	201
Total	108	199	307
Positive Agreement (95% CI): 104/108 96.30% (90.86% - 98.55%)			
Negative Agreement (95% CI): 197/199 98.99% (96.41% - 99.72%)			
Total coincidence rate (95% CI): (104+197)/307 98.05% (95.80% - 99.10%)			

The data below is for understanding information refer to RT-PCR cycle threshold (Ct):

The performance of the SARS-CoV-2 Antigen Test Kit with positive results stratified by the RT-PCR method Ct counts was assessed to more understand the correlation of assay performance to the RT-PCR Ct value, estimating the viral amount present in the clinical sample. As shown in the following table, the positive agreement of the SARS-CoV-2 Antigen Test Kit is higher with samples of a Ct count <33.

Performance against the RT-PCR Method – by Ct Counts

SARS-CoV-2 Antigen Test Kit	RT-PCR Method (Ct)	
	POS (Ct < 33)	POS (Ct ≥ 33)
Positive	95	15
Negative	0	6
Total	95	21
Positive Agreement (95% CI) 100.0 (96.0, 100.0) 71.4 (43.5, 87.4)		

Endogenous Interfering Substances

The following interfering substances, that may be introduced into the nasal cavity or nasopharynx, were evaluated with the SARS-CoV-2 Antigen Test Kit at the concentrations listed below and were no effect on the performance of the test kit.

S.N.	Medicine Name	Concentration Tested	S.N.	Medicine Name	Concentration Tested
1	α-Interferon	10 million U/mL	13	Levofloxacin	500 mg/mL
2	Zanamivir	50 mg/mL	14	Azithromycin	1 g/mL
3	Ribavirin	2 g/mL	15	Ceftriaxone	2 g/mL

4	Oseltamivir	200 mg/mL	16	Meropenem	2 g/mL
5	Paramive	1 g/mL	17	Tobramycin	1 g/mL
6	Lopinavir	1 g/mL	18	Phenylephrine	50 mg/mL
7	Ritonavir	250 mg/mL	19	Oxymetazoline	0.5 mg/mL
8	Abidor	1 g/mL	20	Beclomethason	2 mg/mL
9	Dexamethasone	20 mg/mL	21	Flunisolide	5 mg/mL
10	Triamcinolone acetate	100 mg/mL	22	Budesonide	2 mg/mL
11	Mometason	1 mg/mL	23	Fluticasone	10 mg/mL
12	Histamine hydrochloride	5 mg/mL	24	Sodium chloride (with preservatives)	10 µg/mL (Benzalkonium chloride 50 µg/mL)

LIMITATIONS

1. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly.
2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
3. Test results must be evaluated in conjunction with other clinical data available to the physician.
4. This test cannot distinguish between asymptomatic carriers and infected persons of the SARS-CoV-2.
5. A false-negative result may be obtained if the concentration of the viral antigen in the specimen (swab/saliva) is below the sensitivity.
6. Negative results should be treated as presumptive and confirmed with an approved molecular assay.
7. Clinical performance was evaluated with frozen samples, and performance may be different with fresh samples.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on SARS-CoV-2 Antigen Test Kit are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2016.


Company Information (Specimen Collection Swab)

 Hunan Runmei Gene Technology Co., Ltd.

Tel: +86-731-89919680

Registered and manufacturing address: Room 401-1, Building 3, Shanhe Medical and Health Industrial Park, No. 1048, Zhongqing Road, Shaping Street, Kaifu District, Changsha, Hunan Province, P. R. China.


Post code: 410153



Lotus NL B.V.

Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

Company Information (SARS-CoV-2 Antigen Test Cassette)

 Shenzhen Uni-medica Technology Co. Ltd

Tel: 86-755-86502782 **Fax:** 86-755-86936803

Registered and manufacturing address: Room 202, Block 6th, Liuxian Culture Park, Xili Town, Nanshan District, Shenzhen, Guangdong Province, P. R. China.

Post code: 518055

Website: www.uni-medica.com



CMC MEDICAL DEVICES & DRUGS S.L.

Address: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Manufactured By		CE Mark
	Authorized Representative		Catalog Number
	In Vitro Diagnostic Medical Device		Potential Biological Hazards After Use
	Batch Code		Do Not Reuse
	Expiration Date in Year-Month-Day Format		Date of Manufacture
	Temperature Limitation		Consult instructions for use
	caution		Keep away from sunlight

EC REP CERTIFICATE



CMC MEDICAL DEVICES & DRUGS SL NO. CMC/CE/2021/07042021.1

CONFIRMED THAT CMC MEDICAL DEVICES & DRUGS S.L. Is the European Authorized
Representative of

Shenzhen Uni-medica Technology Co. Ltd
Room 202,Block 6th,LiuXian Culture Park, XiLi Town,
NanShan District, ShenZhen, Guangdong Province, P.R.China

The certificate remains valid until the expiration agreement of EC REP, manufacturing conditions, the quality system or relevant legislation are changed. The validity is conditioned by positive results of periodic surveillance audits.

The product liability rests with the manufacturer in accordance with applicable directive and standard, after fulfilling of the relevant EU legislation requirements, the manufacturer shall affix relevant CE marking to all above mentioned models of the medical device.

Complies with the applicable essential requirements of the council directive 98/79/EEC on in vitro diagnostics as amended.

The products in Annex I was registered in Spanish MOH with number **RPS/659/2021**



Issued on: 07/04/2021

Valid until: 08/06/2022


Authorized Signatory
CMC Medical Devices & Drugs SL


EC REP CERTIFICATE



ANNEX I Medical Device Products



SARS-CoV-2 Antigen Test Kit (Colloidal Gold)





EC Declaration of Conformity

according to the Directive 98/79/EC
(applicable to IVD Devices of NOT Annex II and NOT self-test)



Manufacturer: Shenzhen Uni-medica Technology Co. Ltd

Address: Room 202, Block 6th, Liuxian Culture Park, Xili Town, Nanshan District, Shenzhen, Guangdong Province, P.R.China.

EC Representative: CMC MEDICAL DEVICES & DRUGS S.L.
C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain

We, the manufacturer, declare under our sole responsibility that

Product Name SARS-CoV-2 Antigen Test Kit (Colloidal Gold)

Type 20 Tests/kit, 1 Test/kit

Classification Common/Others IVD
(Devices of NOT Annex II and NOT self-test)

The medical device is in conformity with the relevant provisions and requirements of Directive 98/79/EC of the European Parliament and of the Council on In Vitro Diagnostic Medical Devices.

The products comply with the essential requirements in accordance with Annex I of the In Vitro Devices Directive 98/79/EEC.

Applied harmonised standards, national standards or other normative documents

EN ISO 15223-1:2016
EN 1041:2008
EN ISO 14971:2012
EN 13612:2002
EN 23640:2015

Conformity assessment procedure **Module A (EC Declaration of Conformity) (Annex III, except point 6)**

Signed on: 25th March, 2021 Place: Shenzhen, Guangdong Province, China

Signature (on behalf of the manufacturer):

Name of authorized signatory: 郭永超 GUO YONGCHAO

Position held in the company: General Manager

Company Seal/Stamp:



1 - Product and Company Information

Product Name	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)
Company	Shenzhen Uni-medica Technology Co. Ltd Room 202, Block 6th, Liuxian Culture Park, Xili Town, Nanshan District, Shenzhen, Guangdong Province, P. R. China
Telephone	+86-0755-86505501
General use	Reagent is applied for human diagnostic uses.

2 - Hazards Identification**Classification and Labeling according to Regulation (EC) No 1272/2008 (CLP)**

Hazard category	Not a hazardous substance or mixture according to Regulation (EC)No. 1272/2008		
Pictograms and signal word	None	Hazard statements	None
Precautionary statements	None	Supplemental Hazard Informations (EU)	None
Other information	None		

3 - Composition/Information on Ingredients

Substance/Mixture : Mixture

Components	CAS#	EC#
Test card	No data available	No data available
sample buffer.	7647-14-5	No data available
swab	No data available	No data available

4 - First Aid Measures

After Inhalation

If inhaled, remove to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms, consult a physician.

After skin contact

In case of contact, immediately wash skin with soap and copious amounts of water. If experiencing skin symptoms, consult a physician

After eye contact

In case of contact, immediately flush eyes with copious amounts of water for at least 5 minutes. If experiencing eye symptoms consult a physician

After Ingestion

If swallowed, vomiting. Give large quantities of water and rinse mouth. If feel unwell, consult a physician

5 - Fire Fighting Measures

Extinguishing Media: This product is not combustible. Use any Media appropriate for the surrounding fire.

6 - Accidental Release Measures

Wear appropriate gloves/protective clothing/eye protection/breathing apparatus

Soak up the leakage with inert absorbent material and recover into the suitable and closed containers for disposal.

7 - Handling and Storage

Handling

Directions for Safe Handling: Do not breathe vapor or drink. Avoid contact with eyes, skin, and clothing.

Avoid prolonged or repeated exposure.

Keep away from sources of heat.

Storage

Stable at room temperature and expires in 12 months.

8 - Exposure Controls / Personal Protection

Engineering controls	In general, dilution ventilation is a satisfactory health hazard control for this substance however, if conditions of use create discomfort to the worker, a local exhaust system should be considered Maintain eye wash fountain and quick-drench facilities in work area.
Inhalation protection	Suitable respiratory protective device recommended
Eye protection	Use safety goggles
Skin protection	Wear general protective gloves and clothings
Other protection	Do not eat, drink or smoke in work area

9 - Physical and Chemical Properties

Appearance and odor	Liquid	PH	Not available
Freezing point(°C)	Not available	Boiling point(°C)	Not available
Density (water=1)	No data available	Relative vapour density (air=1)	No data available
Vapour pressure(kpa)	No data available	Heat of combustion (KJ/mol)	No data available
Critical temperature(°C)	No data available	Critical pressure (Mpa)	No data available
Octanol/water partition coefficient as log pow	No data available	Flash point (°C)	No data available
auto-ignition temperature(°C)	No data available	Solubility	Miscible with water
Upper explosive limits%(V/V)	No data available	Lower explosive limits %(V/V)	No data available
Other properties	No data available		

10 - Stability and Reactivity

Stability

Stable under ordinary conditions of use and storage

Incompatible Materials

Strong acids, strong bases, strong oxidizers

Conditions to Avoid	Heat or flame
Hazardous polymerization	Will not occur
Hazardous decomposition	No data available

11 - Toxicological Information

Acute toxicity	No data available		
Skin irritation/corrosion	No data available	Respiratory or skin sensitization	No data available
Eye damage/irritation	No data available	Reproductive cell mutagenicity	No data available
Reproductive toxicity	No data available	Aspiration hazard	No data available
STOT-single exposure	No data available	STOT-repeated exposure	No data available
Carcinogenicity	This material is not listed in IARC (International Agency for Research on Cancer) Category		
Health hazards	Inhalation: May cause mild irritation. May be harmful if inhaled in excessive quantities Ingestion: May be harmful if swallowed Skin Contact: May cause mild irritation or allergic skin reaction Eye Contact: May cause mild irritation		
Other information	No data available		

12 - Ecological Information

Ecological toxicity	This material is not expected toxic to aquatic life
Persistence and degradability	No data available
Bioaccumulation	No data available
Mobility in soil	No data available
PBT and vPvB assessment	No data available
Other information	No data available

13 - Disposal Considerations

Reagent is not a regulated hazardous material.

It is water soluble and may be disposed of as aqueous waste in accordance with state and local regulatory requirements.

14 - Transport Information

Regulations	IATA DGR	IMDG CODE
UN No.	Not regulated	Not regulated
Proper shipping Name	Not regulated	Not regulated
Hazard Class	Not regulated	Not regulated
Packing Group	Not regulated	Not regulated
Packing Method	Not regulated	Not regulated
Environmental Hazards	The substance does not belong to environmentally hazardous substance	
Note	Non-hazardous for air transport	

15 - Regulatory Information

Chinese regulations

Regulations on the Safety Administration of Dangerous Chemicals(2011)

This material was not listed in General rule for classification and hazard communication of chemicals(GB 13690-2009)

EU regulations

Commission Regulation(EC) No. 1907/2006(REACH)and the amendments

Commission Regulation(EC)No, 1272/2008(CLP) and the amendments

Waste Framework Directive 2008/98/EC and the amendments

16 - Other Information

WARRANTY

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

This substance is not dangerous under current provisions of the code of International Carriage of Dangerous Good by Road (ADR) and by Rail (RID),of the International Maritime Dangerous Goods Code (IMDG),and of the International Air Transport Association(IATA DGR 62ND edition) regulations.

Certificate of Analysis

Product Name	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)				
Specs	20 tests/kit	Lot	20210801	Manufacture Date	2021.08.03
Extraction	罗梅莹	Analyst	周冬云	Test Date	2021.08.03
Standard	Product Standard of SARS-CoV-2 Antigen Test Kit (Colloidal Gold)				
QC Result					
Test Item	Standard		Result		Conclusion
Appearance Characteristics	Appearance of kit should be in good condition. Information of all labels should be clear. Contents of Kits should be correct.		Appearance: passed. Information: passed. Contents of Kits: passed.		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Physical examination-Liquid velocity	≥5mm/min		6.4mm/min		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Accuracy-positive coincidence	100.0% while testing positive controls.		100%		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Accuracy-negative coincidence	100.0% while testing negative controls.		100%		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Sensitivity	Test sensitivity control L3, L2 and L1 (3 replicates for each control), the results of testing L1 and L2 should be positive and that of L3 should be negative.		L1: Positive L2: Positive L3: Negative		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Precision	Test control J1 and J2 (10 replicates for each control), the results of testing should be positive.		J1: Positive J2: Positive		<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed
Analyst	周冬云		Date	2021.08.03	
QC Review	张莹		Date	2021.08.03	
Conclusion	<input checked="" type="checkbox"/> passed <input type="checkbox"/> failed				
	QA Review	马强	Date	2021.8.3	
	Certified by	陈文乐	Date	2021.8.3	